


Product Data Sheet

Product name	Medlogics Anti - Decubitus Seat Cushion
REF	AC100
General Information	The Viscoelastic seat cushion offers residents a high level of comfort. It adapts to the individual body shape and provides optimum pressure relief. The Viscoelastic cushion is used for prevention of Decubitus.
Specification	3 layers, with wash label, with insert card, up to 100 kg
Color	Blue
Size	Length 45 cm
	Width 45 cm
	Height 9 cm
Material	Top layer - 3 cm Visco-elastic foam (85 kg/m ³)
	Middle layer - 3 cm soft support foam (28 kg/m ³)
	Bottom layer - 3 cm solid support foam (30 kg/m ³)

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

