

# Product Data Sheet

<b>Product name</b>	Medlogics Bed Pads With Wings
<b>REF</b>	ADL610
<b>General Information</b>	The incontinence bed pads offer optimum protection as well as comfort. The reusable bed pads are made up of four functional layers. The top layer transports unwanted moisture directly to the layers below and dries up in a matter of minutes so that the user will not feel any moisture at all.
<b>Specifications</b>	4 layers, with wings
<b>Color</b>	Green (Top)
	Blue (Bottom)
<b>Size</b>	90 x 85 cm (wings 45 x 85 cm)
<b>Materials</b>	1 <sup>st</sup> layer - 100% brushed polyester (85 g/m <sup>2</sup> )
	2 <sup>nd</sup> layer - 70% polyester / 30% rayon (340 g/m <sup>2</sup> )
	3 <sup>rd</sup> layer - PU laminating (30 g/m <sup>2</sup> )
	4 <sup>th</sup> layer - 100% brushed polyester (85 g/m <sup>2</sup> )
	Wings - 50% cotton / 50% polyester (100 g/m <sup>2</sup> )

## Product photo



### Headquarters:

Medlogics.com LTD  
Lyuben Karavelov 6  
9000 Varna, Bulgaria

**E** info@medlogics.com  
**W** www.medlogics.com

