

Product Data Sheet

Product name	Medlogics Disposable Bed Underpads
REF	ADL700
General Information	The disposable bed underpads are highly absorbent and made of 7 - layers. They are a perfect way to save time and energy on washing. A completely waterproof back ensures no liquid can pass through the underpad. Its soft layers are gentle to the skin. They ensure fast absorption and dispersion across the entire pad.
Color	White (Top)
	Green (Bottom)
Size	60 x 90 cm
Main Fabric	Facing - 100% cotton, 100% polyester, cotton/polyester, bamboo
	Soaker - Polyester/Rayon blend non-woven
	Barrier - PU or PVC
	Backing - 100% Polyester

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

