

Product Data Sheet

Product name	Medlogics Disposable bibs
REF	DP100
General Information	Our disposable dual-layer bibs offer optimum protection of clothing. The practical and large bag of the patient bibs, made of especially strong materials, catches all spilled food. The plastic top layer prevents fluids and food from leaking onto patients' clothes, while the soft paper layer underneath ensures the bib is comfortable to wear.
Specification	Practical pocket to catch spills
	No cleaning costs
	Can be tied around the neck
Size	37 x 66 cm
Material	Soft paper (Top layer)
	Polyethylene plastic underlayer (Bottom layer)

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

