

Product Data Sheet

page 1/2

Product name	Medlogics Gelbodies
Protector REF / Size / Dimensions	Heel & Elbow GB010 / S / 20-25 cm GB011 / M / 25-30 cm GB012 / L / 30-36 cm GB013 / XL / 36-41 cm 
	Shin GB021 / M / 25-33 cm GB041 / L / 33-41 cm 
	Heel & Ankle GB030 / S / 20-25 cm GB031 / M / 25-30 cm GB032 / L / 30-36 cm 
	Knee GB040 / S / 23-28 cm GB041 / M / 28-33 cm GB042 / L / 33-41 cm 
	Hip GB050 / S / 71-79 cm GB051 / M / 79-85 cm GB052 / L / 85-93 cm 

Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com



Product Data Sheet

Product name	Medlogics Gelbodies	page 2/2
General Information	Gelbodies provide protection for areas at risk from pressure sores and nerve damage. GelBodies are suitable for Long Term Care and Rehabilitation patients to assist in the prevention of pressure sores. They are used in the hospital, nursing homes and homecare environment worldwide. The Coolband material feels great, is skin compatible, breathable and can be machine washed up to 40°C. The Gelpads can be cleaned with a damp cloth using a standard cleaning agent. The Gelpads should not be machine washed. Gelbodies can be worn for up to 12 hours in any 24 hour period.	
Packaging:	2 x Silicone gel pads 2 x Textile coolbands 1 x A4 flyer (Applies for each article number)	
Material Gelpad	PU Gel; Outside - TPU film	
Material Textile	85% Cotton, 10% Lycra, 5% Polyester	

Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

