

Product Data Sheet

| | |
|------------------------------|--|
| Product name | Medlogics Lift and Recliner Chair |
| REF | RC100 (Blue) RC110 (Red) RV120 (Yellow) |
| General Information | <p>The chair can be raised and leaned back. The lifting function of the chair is specially designed for older people or people with knee problems. The chairlift in the top position is helpful so that people can get up from the chair. The reclining function is mainly for relaxing, such as watching TV, sleeping and so on. The fabric cover is available in blue, red and yellow. Hardly inflammable. The dead weight is 30 kg and the maximum load capacity is 130 kg. The adjustment of the back and leg rest is made by the integrated hand control. 1 engine. Dimensions: seat width 50 cm, seat depth 50 cm, seat height 49 cm, armrest height 61 cm, backrest width 54 cm, backrest height 65 cm, total height 103 cm, total depth 78 cm.</p> |
| Specifications | Hardly inflammable |
| | Adjustable back and leg rest |
| | Manual operation |
| | 1 engine |
| Material | Polyester |
| Maximum load capacity | 130 kg |

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

