

Product Data Sheet

Product name	Medlogics Viscoelastic Mattress Basic 50 kg/m ³
REF	VM100
General Information	Decubitus ulcers, more commonly known as bedsores, are a common and potentially serious problem for bedridden residents. This mattress is a good choice if the residents are exposed to risk of bedsores and need particular protection and treatment against pressure from lying down.
Specification	2 layers, with wash label, with pamphlet
Color	Blue
Size	Length 200 cm
	Width 90 cm
	Height 12 cm
Material	Outer Layer - 100% polyester non-woven (130 kg/m ³)
	Inner Layer - Memory foam (50 kg/m ³) Polyurethane PU foam (30 kg/m ³)

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com



Product Data Sheet

Product name	Medlogics Viscoelastic Mattress Premium 85 kg/m ³
REF	VM200
General Information	Decubitus ulcers, more commonly known as bedsores, are a common and potentially serious problem for bedridden residents. This mattress is a good choice if the residents are exposed to risk of bedsores and need particular protection and treatment against pressure from lying down.
Specification	2 layers, with wash label, with pamphlet
Color	Blue
Size	Length 200 cm
	Width 90 cm
	Height 12 cm
Material	Outer Layer - 100% polyester non-woven (130 kg/m ³)
	Inner Layer - Memory foam (85 kg/m ³) Polyurethane PU foam (30 kg/m ³)

Product photo



Headquarters:

Medlogics.com LTD
Lyuben Karavelov 6
9000 Varna, Bulgaria

E info@medlogics.com
W www.medlogics.com

